# **BRUXOGARD**<sup>™</sup>

### Protects your teeth and jaws from night time bruxing

## HARD VERSION

## **BRUX**OGARD

WWW.MYORESEARCH.COM WORLDWIDE PATENTS, FOR MORE INFORMATION VISIT MYORESEARCH.COM

# **BRUX**OGARD<sup>™</sup>

The **BRUX**OGARD<sup>™</sup> is fitted over the teeth and worn at night while sleeping. It features a unique flat base specifically designed to protect your teeth and jaws from the effects of night tooth grinding (bruxing). The **BRUX**OGARD<sup>™</sup> is best moulded by your dentist. If you need further moulding, follow the below instructions.

#### **BRUX**OGARD<sup>™</sup> HARD MOULDING INSTRUCTIONS READ ALL INSTRUCTIONS BEFORE BEGINNING THE MOULDING.

- 1 Place **BRUX**OGARD<sup>™</sup> (tongue tag down) into bowl of freshly boiled water (100°C / 212°F) for **90 SECONDS.**
- 2 Lift **BRUX**OGARD<sup>™</sup> gently from water using spoon.
- 3 Cool under tap water for ONE SECOND ONLY bringing surface temperature to comfortable level. Quickly move onto next step.
- 4 Watching in mirror, line up centre of **BRUX**OGARD<sup>™</sup> (tongue tag up) with centre line of upper teeth and proceed rapidly with steps 5 to 7.
- 5 Fit **BRUX**OGARD<sup>™</sup> into upper teeth and press firmly into front teeth and molars.
- 6 Place tip of tongue on tongue tag and close mouth while biting down hard on **BRUX**OGARD<sup>™</sup>.

7 Suck in strongly and press edges of **BRUX**OGARD<sup>™</sup> into teeth and gum line through the lips and cheeks (continue for 20 seconds). Remove from mouth and place in a bowl of cold water for 1 minute. Replace into mouth and test for a good firm fit. Repeat steps 1 to 7 if the fit is not accurate (reheat water to 100°C).

To clean, rinse under warm water. Use toothpaste to clean if desired. Heavy bruxers may wear the **BRUX**OGARD<sup>™</sup> out quickly. Replace if splitting and deterioration occurs.



by Myofunctional Research Co. <u>Myofunctional Research Co. Australia:</u> 44 Siganto Drive, Helensvale QLD 4212 Australia Tel: +61 7 5573 5999

<u>USA Representative:</u> Myofunctional Research Co. USA 9267 Charles Smith Avenue, Rancho Cucamonga CA 91730 USA PO Box 2781, Rancho Cucamonga CA 91730 USA Tel: + 18 66 550 4696



European Representative: Myofunctional Research Europe B.V. Gompenstraat 21c 5145 RM Waalwijk The Netherlands P.O. Box 718, 5140 AS Waalwijk The Netherlands Tel: + 31 416 651 696





BRUX\_H\_INS\_0818\_ENG\_V1.0.0